

TCL VMD SCRIPT TO ANALYSE TRAJECTORY and Can be UTILIZED as CUSTAMIZED COMMANDS
in VMD

VERSION : analysis_script_V 1.0

Authour : ANJI BABU KAPAKAYLA
IIT KANPUR, INDIA.
(anjibabu480@gmail.com)

Step 1 : SOURCE analysis.tcl in VMD Tk console
Step 2 : Command_name {Arguments}

After Sourcing analysis.tcl script in VMD console, type following
my_commands show : List out all available commands from this script.
--help command_name : Shows the details of command ,with examples.

=====

TIP OF THE DAY :

Keep this analysis.tcl script in a adirectory called vmd_tcl_scripts/
Mention its path in vmd startup file ~/.vmrc

eg. source /home/anjibabu/vmd_tcl_scripts/analysis.tcl

That's it, here after when ever you open VMD automatically analysis.tcl
will be excecuted.
So, U can start using these commands as default commands.

=====

Following analysis can be done by using this script :

Alignment
RMSD
Distance btw two atoms
Measure Angle
Measure Dihedral
Contact map
Print Interactions within given cutoff
Print No. of waters within given cutoff
Print No. of HBonds within given cutoff
Print residue names and resid for given sel
Molecule details (No of frames, No of atoms, No of waters, Box size)
On pbc_box and off Pbc_box and prints box size
save_pdb
save_image : Renders image using Tachyon
save_movie : creates a movie for given frames (gif file) other formats will
be updated soon.
save_coordinates
--help : which gives details of command
my_commands show : shows availilable commands

=====

Future updates:
save_view : saves the visivilization of current vmd
Ramachandran plot
HB plot
PBC Wrapping
delete_frame : Deletes the frames for given intial and final frames

#=====

AVAILABLE COMMANDS AND DETAILS :
=====

1. Command Name : my_commands show
PURPOSE : Shows all available commands written in this script
USAGE&example : my_commands show
Here you have to use "show" as fixed argument to see commands list.
2. Command name : --help
PURPOSE : Prints the details of any mentioned command
USAGE : --help command_name
EXAMPLE : --help distance
EXAMPLE 2 : --help angle
Prints the details of commands distance and angle.
3. Command Name : save_pdb
PURPOSE : Writes pdb file for given options
USAGE : save_pdb {atomselection} {start_frame} {end_frame} {stride}
{molid} {filename}
EXAMPLE : save_pdb "protein" 5 100 5 top file.pdb
It will write file.pdb for "protein" from frame 5 to 100 by skipping every 5 frames.
4. Command Name : align
PURPOSE : To align two molecule for given molids
USAGE : align {molid1} {molid2}
EXAMPLE : align top 1
5. Command Name : rmsd
PURPOSE : Measures Avg.RMSD & std for given selections
USAGE : rmsd "sel1" "molid1" "sel2" "molid2"
EXAMPLE : rmsd "backbone" 0 "backbone" 1
OUTPUT : Generates data into a file "RMSD.dat".
(IF both the molids are same , RMSD will be calculated by taking zeroth frame as reference)
6. Command Name : distance
PURPOSE : Measures Avg Distance & std for any given 2 atoms
USAGE : distance sel1 sel2
EXAMPLE : distance "serial 3418" "serial 3415"
OUTPUT : Generates data into a file "DISTANCE.dat"
7. Command Name : angle
PURPOSE : Measures Avg. angle & std for any given 3 atoms
USAGE : angle "sel1" "sel2" "sel3"
EXAMPLE : angle "serial 3418" "serial 3415" "serial 3395"
OUTPUT : Generates data into file " ANGLE.dat ".
8. Command Name : dihedral
PURPOSE : Measures Avg Dihedral angle & std for any given 4 atoms
USAGE : dihedral sel1 sel2 sel3 sel4
EXAMPLE : dihedral "serial 3418" "serial 3415" "serial 3412" "serial 3395"
EXAMPLE : dihedral "index 3417" "index 3414" "index 3411" "index 3394"
OUTPUT : Generates data into a file " DIHEDRAL.dat ".
9. Command Name : show_residues
PURPOSE : Prints RESID and corresponding RESNAME for given selection and for given frames

```

USAGE      : show_residues sel molid <initial frame> <final frame>
EXAMPLE    : show_residues "resid 170 to 180" top 0 5
EXAMPLE1   : show_residues "(not resname WAT) and within 3 of resid 235" 0
5 10
EXAMPLE2   : show_residues "protein and hydrophobic " top 1 0
EXAMPLE3   : show_residues "protein and basic" top
" In Example3 frame numbers are optional ,
  if you omits frame numbers it prints all basic residues from zeroth frame"
OUTPUT     : Generates data into a file "RESNAMES.dat ".

```

```

10.Command Name : details
PURPOSE        : Prints required molecular details for given molid
USAGE          : details molid
EXAMPLE        : details 1
EXAMPLE1       : details top

```

```

11.Command Name : count_waters
PURPOSE         : Prints Avg. No of waters of entire trajectory or for selcted
frames
                which are present around given selection
USAGE          : count_waters "sel" molid start_frame end_frame
EXAMPLE        : count_waters "within 5 of resid 235" 0 5 10
OUTPUT         : Generates data into a file "NO-OF-WATERS.dat".

```

```

12.Command Name : pbc_box
PURPOSE         : Sets BOX ON or OFF or shows the size of PBC BOX
USAGE          : pbc_box {choice}
ARGUMENTS      : on, off, size
EXAMPLE        : pbc_box on
( It will ON the pbc box for current top molecule)
EXAMPLE1       : pbc_box off
( It will off the PBC box )
EXAMPLE2       : pbc_box size
( It prints the size of BOX (values) on Tk console)

```

```

13.Command Name : show_interactions
PURPOSE         : Prints Interactions between given two selections within given
DISTANCE cutoff
USAGE          : show_interactions sel1 sel2 cutoff <start frame> <end frame>
EXAMPLE        : show_interactions "protein" "resid 143" 3.0 5 10
( It prints all interactions within cutoff of 3.0 (A) between protein and
resid 143 of all atoms )
OUTPUT FORMAT:indexA residA resnameA(typeA)---indexB residB
resnameB(typeB) : Distance
ARGUMENTS      :
  sel1          : Any given atom selction
  sel2          : Any given atom selction
  cutoff        : distance cutoff
  start frame   : Starting frame , from here start analysing
  end frame     : Final framme , It stops here
OUTPUT         : Stores all the information in "INTERACTIONS.dat" file

```

```

14.Command Name : contact_map
PURPOSE         : Measures all the contacts between given two selections and
                  Prints 1 (one) if contact is below cutoff , 0 (zero)
otherwise
USAGE          : contact_map " sel1" "sel2" cutoff startframe endframe
Example        : contact_map "resid 1 to 50 and name CA" "resid 51 to 100 and
name CA" 6.0 0 10

```

```

OUTPUT      : Stores all data in "CONTACT_MAP.dat" file.
OUTPUT FORMAT : RESID_A  RESID_B 1/0 DISTANCE

PLOTTING DATA : gnuplot >>  splot "CONTACT_MAP.dat" u 1:2:3 w p notitle
                  gnuplot >>  set view 0,0,1
                  gnuplot >>  replot

```

It generates a beautiful CONTACT MAP 3D image.

```

15.Command Name : count_hbonds
PURPOSE        : Prints AVG No. of HBonds b/w given DONOR and ACCEPTOR
selections
                  and given distance_cutoff & angle_cutoff .
USAGE          : count_hbonds "Donor_sel1" "Acceptor_sel2" distance_cutoff
angle_cutoff
EXAMPLE        : count_hbonds "protein and name N" "protein and name O" 3.0
30
OUTPUT         : Stores data in file "HBONDS_COUNT.dat".
(Which contains Frame vs No.of Hbonds ,You can plot this data as Frames v/s
No. of Hbonds)
NOTE          : It assumes that molis as "top" or zero. It prints AVG HBONDS
per FRAME.

```

```

16.Command Name : show_hbonds
PURPOSE        : Prints the all HBonds b/w given DONOR & ACCEPTOR selections
for
                  given distance & angle cutoff.
USAGE          : show_hbonds "D_sel1" "A_sel2" dist_cutoff angle_cutoff
<startframe> <endframe>
EXAMPLE        : show_hbonds "protein and name N" "protein and name O" 3.0 30
5 10
OUTPUT         : Stores all HBONDS data in HBONDS.dat.
OUTPUT FORMAT  : indexD residD resnameD(type)-- indexA resida
resnameA(type)-- Distance
(Prints and Stores data for every Frame )
NOTE          : Always Atomselctons should be DONOR and ACCEPTOR , (not
ACCEPTOR and DONOR)

```

```

17.Command Name : HB_occupancy
PURPOSE        : Prints only the HBONDS occupancy for given selctions for
selected frames
USAGE          : HB_occupancy "D_sel1" "A_sel2" dist_cutoff angle_cutoff
<startframe> <endframe>
EXAMPLE        : HB_occupancy "protein" "resid 147" 3.0 30 0 20
OUTPUT FORMAT  : RESID_D  RESNAME_D --  RESID_A  RESNAME_A    ---OCUPANCY
(%)
EXAMPLE        : 140  LYS -- 147 ARG  --- 80.00 %    (Pints for all selected
frames)
ARGUMENTS      :
SELECTION1     : MUST BE ANY DONOR SELECTION
SELECTION2     : MUST BE ANY ACCEPTOR SELECTION
CUTOFF         : Distance cutoff to measure HBonds
ANGLEE_CUTOFF : Angle cutoff to measure HBonds
START_FRAME    : Starting Frame Number
END_FRAME      : End Frame Number

```

```

18.Command Name : save_image
PURPOSE        : Renders a high quality image using Tachyon for currently
active frame.
UASAGE        : save_image {molid} filename
EXAMPLE       : save_image top my_image

```

(extension of file not required)
 OUTPUT FILES : filename.dat , filename.dat.tga, filename.dat.jpg
 Above example produces my_image.dat , my_image.dat.tga, & my_image.dat.jpg.

19.Command Name : save_movie
 PURPOSE : GENERATES A MOVIE IN GIF FORMAT FOR GIVEN MOLID & FOR
 SELECTED FRAMES
 USAGE : save_movie molid filename <start frame> <end frame>
 (no extension of filename is required)
 EXAMPLE1 : save_movie top my_protein 0 20
 EXAMPLE2 : save_movie 1 drug_protein 20 40
 EXAMPLE3 : save_movie top my_protein
 EXAMPLE4 : save_movie 1 drug_protein 30
 OUTPUT : Example1 Genarates "my_protein.gif" file from zeroth frame
 to 20th frame.
 : Example2 generates "drug_protien.gif" file from 20th frame
 to 40th frame.
 : Example3 generates "my_protein.gif" file from zeroth frame
 to end frame .
 : Example4 generates "drug_protein.gif" file from 30th frame
 to end frame.

NOTE : Other Movie formats will be updated soon.

20.Command Name : save_coordinates
 PURPOSE : Saves coordinates for given arguments
 USAGE : save_coordinates "Sel" {molid} {filename} {filetype}
 {start_frame} {end_frame} {stride}
 EXAMPLE : save_coordinates "protein" top my_protien.mol2 mol2 5 30 5
 ARGUMENTS :
 Atomselection: Any given slection
 Start frame : From which frame to save coordinates
 End frame : To which frame to save coordinates
 Stride : After how many frames you want save (skipping)"
 molid : Which molecule you want save
 Filename : Filename with extension
 EXAMPLE : protein.mol2, protein.xyz, protein.gro, protein.dcd.....etc.
 Filetype : Which format you want to save.

Available Filetypes : ABINIT , bgf, binpos, crd, crdbox, dcd, gro, trr,
 js, lammprj, mol2, namdbin, pdb, pqr, rst7,POSCAR, xbgf, xyz, dtr, m"ae, dms,
 hoomd .

```

=====
#
#                               WRITTEN BY ANJI BABU KAPAKAYALA
#                               INDIAN INSTITUTE OF TECHNOLOGY,KANPUR.
#                               (anjibabu480@gmail.com)
#
#=====

```

" yad bhavam tad bhavati " means "you become whatever you think"

----from vedas

